


# HIPP Endoskop Service GmbH

>> your requirements are our aims.



*Technische Produktinformation LED-Lichtquellen*  
*Technical Product Information LED Light Sources*






- LED Technologie
- Vergleichbar mit einer 300 Watt Xenon-Lichtquelle
- LED technology
- comparable to a 300 watt xenon light source

Technical Data

Lichtleiteranschluss  
 Light guide connection:




Lampe Lamp	Daylight quality, 6000 K (+/-500 K) Lifetime: ≤30.000 hr
Funktionen Features	automatic light intensity control
Energieversorgung Power supply	Voltage range: 100-240 V, power consumption 130 W Device protection: 2 x T6,3 A Frequency:50-60 Hz
Video IN	1x Video
Standards	Protection class (acc. to IEC 60601-1: I Safety class: BF IP code: 20 Standard: ICE 60601-1-2 Device classification acc. MDD 93/42/EEC: I
Abmessungen Dimension	WxHxD: 355mm x 120mm x 340mm wight: ca. 10 kg
Lieferumfang Delivery contains	Light source, mains cable, user manual, BNC cable



- LED Technologie
- Vergleichbar mit einer 180 Watt Xenon-Lichtquelle
- LED technology
- comparable to a 180 watt xenon light source

Technical Data

Lichtleiteranschluss  
 Light guide connection:



LED	LED CRI 70, Daylight quality, approx. 6500 K Lifetime: > 20.000 hr
Funktionen Features	automatic light intensity control
Energieversorgung Power supply	Voltage range: 100-240 V, power consumption 130 W Device protection: 2 x T6,15 A Frequency:50-60 Hz
Anschlüsse Connections	Video IN: 1x Video, BNC Service (RJ45)
Standards	Protection class (acc. to IEC 60601-1: I Safety class: BF IP code: 20 Standard: ICE 60601-1-2 Device classification acc. MDD 93/42/EEC: I
Abmessungen Dimension	WxHxD: 360mm x 124mm x 335mm wight: ca. 5 kg
Lieferumfang Delivery contains	Light source, mains cable, user manual, BNC cable



#### Product Profile

- |  |  |
|--|--|
| <ul style="list-style-type: none"><li>• 5 einstellbare Helligkeitsstufen</li><li>• Batterie mit langer Lebensdauer</li><li>• höchste Kompatibilität durch ACMI-Anschluss garantiert</li><li>• kein Lampenwechsel dank LED-Technik</li><li>• 37 Minuten Betriebsdauer bei max. Beleuchtungsstärke</li><li>• wasserdicht, komplett einlegbar</li><li>• inkl. Magnetischem USB ladechip</li></ul> | <ul style="list-style-type: none"><li>• 5 adjustable brightness levels</li><li>• long life time battery</li><li>• ultimate compatibility assured by ACMI connection</li><li>• no exchange of lamps necessary due to LED technique</li><li>• 37 minutes non-stop operation at maximum power</li><li>• waterproof, completely submersible</li><li>• incl.USB-Magnetic chip</li></ul> |
|--|--|

#### Technical Data

Model / model / Article Nr.	MobiLED / HT-LED3-ZU
charge time	2 h
Watt / watt	10 Watt
Größe / length	63mm x 21mm
Gewicht / weight	51g inkl. Batterie/51g including battery
Helligkeit / Brightness	Max 1000 Lumen / max lumens
Helligkeitsstufen / brightness level	Turbo / High / Medium / Low / Moon
Lieferumfang / Delivery contains	Lichtquelle / Ladestation / USB-Ladekabel / Gebrauchsanweisung / Light source / Charging station / USB charging cable / Instructions for use

