


HIPP Endoskop Service GmbH


>> *your requirements are our aims.*



Technische Produktinformation Hysteroskope
Technical Product Information Hysteroscopes



Hysteroskope Hysteroscopes



Hysteroskope – Storz-kompatibel
Hysteroscopes – Storz-compatible

autoklavierbar (134°C / 273°F)
autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
2,9	0°	302	x	30-0705-00
2,9	30°↑	302	x	30-0707-00
4,0	0°	302	x	30-0731-00
4,0	12°↑	302	x	30-0064-00
4,0	30°↑	302	x	30-0775-00
4,0	70°↑	302	x	30-0066-00

