

HIPP Endoskop Service GmbH


>> Your requirements are our aims




Technische Produktinformationen Arthroskope
Technical product information Arthroscopes



*Made in
Germany*






Arthroskope Arthroscopes

Arthroskope- Storz- kompatibel Arthroscopes - Storz- compatible				autoklavierbar (134° / 273°F) autoclavable (134°C / 273°F)
Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article no.
1,9	0°	60	x	30-0546-00
1,9	30°	60	x	30-0433-00
1,9	0°	102	x	30-0256-00
1,9	30°	102	x	30-0255-00
2,7	0°	110	x	30-0793-00
2,7	30°	110	x	30-0752-00
2,7	70°	110	x	30-0003-00
2,7	0°	141	x	30-0265-00
2,7	30°	141	x	30-0784-00
2,7	0°	175	x	30-0708-00
2,7	30°	175	x	30-0709-00
2,7	70°	175	x	30-0710-00
2,7	0°	187,5	x	30-0742-00
2,7	30°	187,5	x	30-0743-00
2,7	70°	187,5	x	30-0744-00
4,0	0°	175	x	30-0693-00
4,0	30°	175	x	30-0694-00
4,0	45°	175	x	30-0695-00
4,0	70°	175	x	30-0696-00





Arthroskope - Dyonics- Kompatibel Arthroscopes - Dyonics- compatible				autoklavierbar (134° / 273°F) autoclavable (134°C / 273°F)
Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article no.
4,0	30°	160	x	30-0725-00
4,0	70°	160	x	30-0765-00

Arthroskope - Stryker- Kompatibel Arthroscopes - Stryker- compatible				autoklavierbar (134° / 273°F) autoclavable (134°C / 273°F)
Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article no.
4	30°	141	x	30-0324-00
4	70°	141	x	30-0344-00

Arthroskope - Wolf- kompatibel Arthroscopes - Wolf- compatible				autoklavierbar (134° / 273°F) autoclavable (134°C / 273°F)
Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article no.
4	25°	170	x	30-0755-00
4	30°	175	x	30-0429-00
4	70°	177	x	30-0668-00

