


HIPP Endoskop Service GmbH


>> *your requirements are our aims.*



Technische Produktinformation Arthroscopie
Technical Product Information Arthroscopes




Arthroskope Arthroscopes



Arthroskope – Storz-kompatibel
Arthroscopes – Storz-compatible

autoklavierbar (134°C / 273°F)
autoclavable (134°C / 273°F)


Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
1,9	0°	60	x	30-0546-00
1,9	30°	60	x	30-0433-00
2,7	0°	110	x	30-0793-00
2,7	30°	110	x	30-0752-00
2,7	70°	110	x	30-0003-00
2,7	0°	141	x	30-0265-00
2,7	30°	141	x	30-0784-00
2,7	0°	175	x	30-0708-00
2,7	30°	175	x	30-0709-00
2,7	70°	175	x	30-0710-00
2,7	0°	187,5	x	30-0742-00
2,7	30°	187,5	x	30-0743-00
2,7	70°	187,5	x	30-0744-00
4,0	0°	175	x	30-0693-00
4,0	30°	175	x	30-0694-00
4,0	45°	175	x	30-0695-00
4,0	70°	175	x	30-0696-00



Arthroskope – Dyonics-kompatibel
Arthroscopes – Dyonics-compatible

autoklavierbar (134°C / 273°F)
 autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
4.0	30°	160	x	30-0725-00
4.0	70°	160	x	30-0765-00



Arthroskope – Stryker-kompatibel
Arthroscopes – Stryker-compatible

autoklavierbar (134°C / 273°F)
 autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
4	30°	141	x	30-0324-00
4	70°	141	x	30-0344-00

