


HIPP Endoskop Service GmbH


>> *your requirements are our aims.*



Technische Produktinformation Cystoskope
Technical Product Information Cystoscopes




Cystoskope Cystoscopes



Cystoskope – Storz-kompatibel
Cystoscopes – Storz-compatible

autoklavierbar (134°C / 273°F)
autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
2,7	0°	187,5	x	30-0269-00
2,7	30°	187,5	x	30-0268-00
2,7	0°	302	x	30-0053-00
2,7	30°	302	x	30-0054-00
2,7	70°	302	x	30-0134-00
2,9	0°	302	x	30-0704-00
2,9	30°	302	x	30-0706-00
4,0	0°	302	x	30-0750-00
4,0	12°	302	x	30-0714-00
4,0	30°	302	x	30-0715-00
4,0	70°	302	x	30-0716-00



30° ↑

Cystoskope – Olympus-kompatibel
Cystoscopes – Olympus-compatible

autoklavierbar (134°C / 273°F)
autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
3,0	30°	282	x	30-0510-00

